



IODONTOSUL
PRODUCTS FOR HEALTH

CATALOG
2024

HIGHLIGHTS







For 27 years, IODONTOSUL has been dedicated to manufacturing high-quality dental products to meet the needs of professionals in the field. With a team of highly qualified experts and a commitment to excellence, we have become leaders in our segment. Our products are carefully designed and manufactured with the best materials available on the market, aiming to guarantee the safety and comfort of patients, as well as the effectiveness of treatments. We understand that dentistry is a constantly evolving science, so we are always up to date with the latest technologies and trends in our sector.

All of our products undergo rigorous quality testing to ensure they meet the highest performance standards.

We are committed to providing exceptional service to our customers, working closely with them to understand their expectations and providing customized solutions to meet those needs.

If you are a dental professional looking for high quality and reliable products, do not hesitate to contact us. We will be happy to help you find the perfect solution to suit you.

Our **greatest asset**
is having **you as a**
customer!



1

FINISHING
AND POLISHING

2

AGENT OF
CONNECTION

3

ACESSÓRIOS

4

8

AGENTS
CONDITIONING

9

PREVENTION

10

RADIOLOGY

11

5

6

7




VARIOUS



CEMENTS



WHITENING



ENDODONTIC

12

13

14




RESTORATOR
TEMPORARY



PIN
FIBERGLASS



REMOVER
FROM STAIN



LUBRICANT

Baking soda

CODE. BAR	DESCRIPTION	UN	I.P.I	NCM
7898936469717	Natural Bottle 100g	-	-	2836.30.00
7898936469724	Natural Bottle 250g	-	-	2836.30.00
7898936469731	Natural Bottle 500g	-	-	2836.30.00
7898936469748	Natural box with 15 sachets 40g	-	-	2836.30.00
7898936469755	Natural box with 24 sachets 40g	-	-	2836.30.00



**Quality in the
removal of tartar
and plaque.**

Registro ANVISA: 80442020074



CODE. BAR	DESCRIPTION	UN	I.P.I	NCM
7898936469618	Bottle mint 100g	—	—	2836.30.00
7898936469625	Bottle mint 250g	—	—	2836.30.00
7898936469632	Bottle mint 500g	—	—	2836.30.00
7898936469649	Mint box with 15 sachets 40g	—	—	2836.30.00
7898936469656	Mint box with 24 sachets 40g	—	—	2836.30.00

CODE. BAR	DESCRIPTION	UN	I.P.I	NCM
7898936469663	Bottle strawberry 100g	—	—	2836.30.00
7898936469670	Bottle strawberry 250g	—	—	2836.30.00
7898936469687	Bottle strawberry 500g	—	—	2836.30.00
7898936469694	Strawberry box with 15 sachets 40g	—	—	2836.30.00
7898936469700	Strawberry box with 24 sachets 40g	—	—	2836.30.00



PAST PRO

Prophylactic Paste

RECOMMENDATION:

The product is indicated for prophylactic cleaning to remove exogenous stains, plaques, alba matter and oral residues.

FEATURES AND ADVANTAGES:

Does not damage the contra-angle;
Flavors: Mint and Tutti-Frutti.

REG. ANVISA: 80442020050

CODE. BAR	DESCRIPTION	UN	I.P.I	NCM
7898936469229	Mint	12	-	3006.40.12
7898936469205	Tutti-Frutti	12	-	3006.40.12



**Optimal abrasiveness
and patient comfort.**

REG. ANVISA: 80442020008

POMITEC

RECOMMENDATION:

Recommended for finishing and polishing tooth enamel, cohesive gold, amalgam and acrylic resin.

FEATURES AND ADVANTAGES:

Extra-fine abrasive powder with a greyish color, low density and weight.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467584	Pumice stones 100g	12	—	2513.10.00
7898936467577	Pumice stones 1000g	1	—	2513.10.00

REG. ANVISA: 80442020076



PAST X

RECOMMENDATION:

The product indicated for prophylactic cleaning in the removal of exogenous stains, plaques, alba matter and oral residues.

FEATURES AND ADVANTAGES:

Does not damage the contra-angle.

CÓD. DE BARRAS	DESCRIÇÃO	UN	I.P.I	NCM
78986987650480	Past X	12	—	3006.40.12

REG. ANVISA: 80442020058

DIAMOND PRÓ



RECOMMENDATION:

The product is recommended for polishing restorative materials in general, especially porcelain, with the help of felt discs.

FEATURES AND ADVANTAGES:

Excellent polishing and shine on materials such as:

- Resin compounds;
- Porcelain;
- Tooth enamel;
- Produced with micronized diamond particles, with extra-fine grain;
- High brightness capacity.

CÓD. DE BARRAS	DESCRIÇÃO	UN	I.P.I	NCM
7898936469373	Diamond Pró	12	-	3006.40.12

REG. ANVISA: 80442029011

SILANO

RECOMMENDATION:

Indicated for the treatment of internal surfaces of porcelain restorations and fiberglass posts (support Glass)

FEATURES AND ADVANTAGES:

- Increase in the adhesiveness of materials;
- Ease of use and less evaporation;
- Easy application, pre-activated solution;
- Greater adherence.



CÓD. DE BARRAS	DESCRIÇÃO	UN	I.P.I	NCM
7898698760145	Silano	12	-	3006.40.12

FINISHING AND POLISHING / BONDING AGENT

REG. ANVISA: 80442020064

BANDA MATRIZ POLIÉSTER

Matrix band tape made of polyester

RECOMMENDATION:

Recommended for finishing and polishing dental enamel, cohesive gold, amalgam and acrylic resin

FEATURES AND ADVANTAGES:

Extra-fine abrasive powder with a greyish color, low density and weight.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468185	Polyester matrix band	50	—	3701.10.21

REG. ANVISA: 80442020034



TIRA DE LIXA POLÉSTER

Abrasive tape made of polyester

RECOMMENDATION:

Strip for finishing and interproximal polishing of composite resin dental restorations.

FEATURES AND ADVANTAGES:

- It has two parts: one thin and one thick, delimited by a neutral center;
- Excellent interproximal polishing;
- Easy handling.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467867	Sandpaper 50 units	10	—	3701.10.21
7898936467881	Sandpaper 100 units	10	—	3701.10.21
7898936467898	Sandpaper 150 units	10	—	3701.10.21

REG. ANVISA: 80442020026

BANDA MATRIZ

Matrix band tape made of stainless steel

RECOMMENDATION:

Stainless steel matrix band strip for reconstruction of dental interproximal contact.

FEATURES AND ADVANTAGES:

- 0.05mm x 5mm x 50cm;
- 0.05mm x 7mm x 50cm;
- Flexible;
- Easy to remove.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468666	0,05mm x 5mm x 50cm	50	—	9018.90.95
7898936468673	0,05mm x 7mm x 50cm	50	—	9018.90.95

TIRA ABRASIVA DE AÇO

Abrasive tape made of steel

REG. ANVISA: 80442020035

RECOMMENDATION:

Strip for finishing and polishing the interproximal surfaces of dental restorations.

FEATURES AND ADVANTAGES:

- Excellent interproximal wear of the tooth structure;
- Composite resin;
- Acrylic resin;
- Amalgam;
- Glass ionomer;
- Easy handling.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467775	Steel Strip 4mm	10	—	3006.40.12
7898936467843	Steel Strip 6mm	10	—	3006.40.12

REFLEX

Dental Mirror

RECOMMENDATION:

It is a device to assist in the professional's visualization and diagnosis.

FEATURES AND ADVANTAGES:


- Developed in AISI 304 stainless steel;
- Autoclavable 121°C;
- Electronically and manually polished;
- Real image, no shadows.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468734	Individual plan	10	-	700991.00
7898936467980	Plan 12 units	10	-	700991.00
7898936468895	Front surface	10	-	700991.00
7898698760602	Individual Plan - Vintage	10	-	700991.00

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760565	Plan 12 units - Vintage	10	-	700991.00
7898698760596	Individual Plan - Rose	10	-	700991.00
7898698760558	Plan 12 units - Rose	10	-	700991.00
7898698760589	Individual Plan - Gold	10	-	700991.00
7898698760541	Plan 12 units - Gold	10	-	700991.00

REG. ANVISA: Dispensa



CABO REFLEX

Dental mirror cable

RECOMMENDATION:

Odontoscope rod extender

FEATURES AND ADVANTAGES:

- Developed in anodized aluminum in 3 colors or 304 stainless steel;
- Autoclavable 121°C;
- Polished bosses;
- Polished electronically and manually.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468826	Chrome-plated aluminium handle	10	-	6804.22.19
7898936468833	Pink aluminium handle	10	-	6804.22.19
7898936467874	Blue aluminium handle	10	-	6804.22.19

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469441	Stainless steel handle	10	-	6804.22.19
7898936469434	Hexagonal Cable	10	-	6804.22.19

REG. ANVISA: 80442029001

CARBON PAPER

Carbon Paper

RECOMMENDATION:

The product is recommended for polishing restorative materials in general, especially porcelain, with the help of felt discs.

FEATURES AND ADVANTAGES:

- Carbon Paper in blue;
- Bicolor Carbon Paper (blue and red);
- Double side;
- Ideal thickness for occlusions.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467539	Carbon Paper Blue	10	—	3006.40.12
7898936467522	Carbon Paper Bicolor	10	—	3006.40.12

PERIOSTONE

Sharpening Stone

RECOMMENDATION:

Stone recommended for sharpening instruments.

FEATURES AND ADVANTAGES:

- Developed in aluminum oxide;
- Varied dimensions;
- Various formats.

REG. ANVISA: Dispensa



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468802	Pink Triangular	10	—	6804.22.19
7898936468819	White gouge	10	—	6804.22.19
7898936468796	Orange guice	10	—	6804.22.19
7898936469021	Gouge pro max orange	10	—	6804.22.19

ANATOMICAL WEDGE



REG. ANVISA: 80442020020

RECOMMENDATION:

They are indicated for use in restorative procedures using matrices, with the mechanical function of adapting the matrix and/or moving adjacent teeth away. The wedges are dehydrated and polished, and therefore do not require cutting and/or adjustments, ensuring the fixation of the matrices and providing restorations free of excess and with a perfect contour.



FEATURES AND ADVANTAGES:

- Used in class II restoration procedures;
- Symmetrical and anatomical geometry;
- Easy handling with the tweezers;
- Colored with non-toxic pigment;
- They do not stain the restoration.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467836	Coloured 100-pcs.	10	—	3006.40.12
7898936467003	Yellow N 1 - 100 pcs.	10	—	3006.40.12
7898936468000	Colorless - 100 pcs.	10	—	3006.40.12

REG. ANVISA: 80442020027

ABRITEX

Mouth Opener

RECOMMENDATION:

It is a device to help open the patient's mouth, facilitating the execution of clinical procedures.

FEATURES AND ADVANTAGES:

- Developed in 100% silicone;
- Autoclavable 121°C;
- Anatomical.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467416	Abritex Infant	10	—	9018.4999
7898936467409	Abritex Adult	10	—	9018.4999
7898936467676	Abritex Couple	10	—	9018.4999

REG. ANVISA: Dispensa



PONTA INTRAORAL

Addition silicone applicator

RECOMMENDATION:

Odontoscope rod extender.

FEATURES AND ADVANTAGES:

- Developed in color anodized aluminum or 304 stainless steel;
- Blast clavable 121°C;
 - Polished projections;
 - Electronically and manually polished.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467010	Number 1 - containing 12 pcs.	10	—	9018.4999
7898936467041	Number 1 - containing 5 pcs.	10	—	9018.4999
7898936467058	Number 1 - containing 10 pcs	10	—	9018.4999
7898936467065	Number 2 - containing 12 pcs.	10	—	9018.4999
7898936467072	Number 2 - containing 50 pcs.	10	—	9018.4999

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467096	Number 3 - containing 12 pcs.	10	—	9018.4999
7898936467386	Number 3 - containing 50 pcs.	10	—	9018.4999
7898936467430	Number 4 - containing 12 pcs.	10	—	9018.4999
7898936467461	Number 4 - containing 50 pcs.	10	—	9018.4999

REG. ANVISA: Dispensa

PLACA DE VIDRO

Glass plate - Polished and ground

RECOMMENDATION:

Recommended for troweling diverse materials.

FEATURES AND ADVANTAGES:

- Available in polished and unpolished - sandblasted;
- Various thicknesses;
- Lapidated.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467614	6mm polished glass plate	10	-	701790.00
7898936467621	6mm ground glass plate	10	-	701790.00
7898936467591	10mm polished glass plate	10	-	701790.00

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467607	10mm ground glass plate	10	-	701790.00
7898936467638	20mm polished glass plate	10	-	701790.00
7898936467645	20mm ground glass plate	10	-	701790.00

MANDRITEC

For sanding strip, contra-angle and straight piece

REG. ANVISA: Dispensa

RECOMMENDATION:

Dental device adapter.

FEATURES AND ADVANTAGES:

- Developed in AISI 304 stainless steel;
- Electronically and manually polished;
- Autoclavable 121°C
- Polished bosses.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468680	CA Contra-angle	10	-	9018.41.00
7898936467492	CA Sanding disc	10	-	9018.41.00
7898936468772	PM Reinforced handpiece	10	-	9018.41.00
7898936468758	PM Handpiece	10	-	9018.41.00
7898936468840	PM Sanding strip	10	-	9018.41.00

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467959	Kit 50	5	—	9018.90.95
7898936467942	Kit 25	5	—	9018.90.95
7898936468864	50 P matrices	10	—	9018.90.95
7898936467423	50 M matrices	10	—	9018.90.95
7898936468871	50 G matrices	10	—	9018.90.95
7898936467553	50 assorted matrices	10	—	9018.90.95
7898936467911	2 Clips + 8 Protectors	10	—	9018.90.95
7898936467935	20 protectors	10	—	9018.90.95

REG. ANVISA: 80442020033



PRO MATRIX

Universal matrices

RECOMMENDATION:

Class II restorations, allowing the restoration of dental anatomy and the creation of contact points.

FEATURES AND ADVANTAGES:

- System composed of pre-formatted sectional matrices and disposables and clamps for fixing;
- Ease of reconstructing the contact point;
- The format of the matrices favors the reestablishment of dental anatomy.



ENDO LINE

Endodontic ruler

RECOMMENDATION:

Used to apply the working length to the instrument by moving the cursor

FEATURES AND ADVANTAGES:

Ruler with central adaptation for the instrument, facilitating length identification.

REG. ANVISA: Dispensa



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760428	Endodontic ruler	50	—	3701.10.21

REG. ANVISA: Dispensa

TAMBOREL

File holder

RECOMMENDATION:

Designed to easily organize files in the procedure and endodontics and their sterilization.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760435	Tamborel	12	—	3701.10.21

REG. ANVISA: Dispensa

PAPEL TOALHA

Paper towel to dry your hands



RECOMMENDATION:

Recommended for drying your hands.

FEATURES AND ADVANTAGES:

- 100% cellulose;
- Two folds
- Color white;
- Format 20cm x 20cm;
- 1,000 sheets.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468352	Paper towels	10	—	4818.20.00

AVENTAL DESCARTÁVEL

Disposable protective apron

RECOMMENDATION:

Recommended for the protection of patients and healthcare professionals.

FEATURES AND ADVANTAGES:

- Weight 40g/m2;
- Made from TNT;
- Elastic cuffs.

REG. ANVISA: 80442020045



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469083	White apron	10	—	6210.10.00

VARIOUS

BABAPRENDI

Napkin clip

RECOMMENDATION:

Device for holding bibs and napkins when serving customers

FEATURES AND ADVANTAGES:

- 40cm long chain;
- Silicone handle, pink and blue colors;
- Not autoclavable.

REG. ANVISA: Dispensa



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467485	Babaprendi	10	8%	9018.4999
7898936467485	Babaprendi	10	8%	9018.4999
7898936467485	Babaprendi	10	8%	9018.4999

VARIOUS

REG. ANVISA: Dispensa



BRUSH DRILL

Drill cleaning brush

RECOMMENDATION:

Recommended for cleaning drill bits

FEATURES AND ADVANTAGES:

- Developed with bristles corrugated brass;
- Does not damage the drill;
- Aluminum handle.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468789	Drill brush	10	-	9603.2900

REG. ANVISA: Dispensa



PONTA MISTURADORA

Tip for mixing and applying materials

RECOMMENDATION:

Recommended for applying molding materials with a homogeneous mixture in 1:1 cartridges.

FEATURES AND ADVANTAGES:

- Double ultraviolet treatment;
- Does not damage the autoclave;
- Non-sterile.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469137	Mixing tip	5	-	9018.49.19

PROPHY BRUSH E CUP

Prophylaxis brush and rubber cup

RECOMMENDATION:

Recommended for dental prophylaxis with the help of PAST PRO, and can also be used for pre-polishing resin with DIAMOND PRO.

FEATURES AND ADVANTAGES:

PROPHY BRUSH: Virgin white bristles;
 PROPHY CUP: virgin white rubber; AC rod.

REG. ANVISA: 80442020047
 80442020084



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467683	Brushy Prophy - Straight Type CA	50	-	9603.21.00
7898936467720	Brushy Prophy - Cup Type	50	-	9603.21.00
7898936467713	Brush Prophy - Conical type CA	50	-	9603.21.00
7898936467034	Cup Prophy - Rubber cup CA	50	-	9603.21.00

VARIOUS

POTE DAPPEN

Plastic Dappen Pot

REG. ANVISA: Dispensa

RECOMMENDATION:

- pot for handling materials.

FEATURES AND ADVANTAGES:

- It has two cavities for manipulation;
- Autoclavable;
- Colors: Blue and Pink.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760527	Blue Plastic Dappen Pot	12	—	701790.00
7898698760534	Pink Plastic Dappen Pot	12	—	701790.00

REG. ANVISA: Dispensa

ESPÁTULA 142

Spatula for Glass Ionomer



RECOMMENDATION:

Product indicated for mixing Ionomers.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760435	Plastic spatula for glass ionomer n 142	10	—	9018.4999

REG. ANVISA: Dispensa

ESPÁTULA ALGINATO

Spatula for handling Alginate



RECOMMENDATION:

- Recommended for mixing alginate, prior to the molding procedure.

FEATURES AND ADVANTAGES:

- Its format allows for great homogenization and hygiene, for having a smooth surface;
- Colors: Blue and Pink.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760510	Pink alginate plastic spatula	10	—	9018.4999
7898698760503	Blue alginate plastic spatula	10	—	9018.4999

REG. ANVISA: Dispensa

PORTA MOLDEIRA

Tray holder

RECOMMENDATION:

- Suitable for storing trays and mobile devices.

FEATURES AND ADVANTAGES:

- Offers greater protection and hygiene for the storage of trays and mobile devices;



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760398	Tray Door	10	—	4911.9900

VARIOUS

REG. ANVISA: 80442020083

CERA ODONTOLÓGICA

Wax

RECOMMENDATION:

- It is indicated for assistance in the creation of dental prostheses, in adjusting trays and adjusting models in articulators.

FEATURES AND ADVANTAGES:

- Composed of paraffin, polyethylene wax and dyes.
- It is soft and has excellent flexibility, with adhesion.
- For dental use.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760664	7 Wax in slides	10	-	340700.20
7898698760657	9 Wax on slides	10	-	340700.20
7898698760640	Wax Articulation A	10	-	340700.20

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760633	Wax Utility U	10	-	340700.20
7898698760626	Ivory Carving Wax	10	-	340700.20
7898698760619	Normal Carving Wax	10	-	340700.20

VARIOUS

REG. ANVISA: Dispensa



POTE PALADON

Silicone Paladon pot

RECOMMENDATION:

- Pot for handling materials.

FEATURES AND ADVANTAGES:

- It has two cavities for manipulation;
- Autoclavable;
- Color: assorted.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760749	Paladon Silicone Pot - Assorted	10	-	701790.00

REG. ANVISA: Dispensa

POTE DAPPEN

Silicone Dappen Pot

RECOMMENDATION:

Pot for handling materials.

FEATURES AND ADVANTAGES:

- It has two cavities for manipulation;
- Autoclavable;
- Color: assorted.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760718	Dappen Silicone P Pot - Assorted	10	-	701790.00
7898698760725	Dappen Silicone Pot G - Assorted	10	-	701790.00

MARCADOR DE INSTRUMENTAL

Instrument Marker

RECOMMENDATION:

Suitable for marking and identifying instruments.

FEATURES AND ADVANTAGES:

- Material: Silicone.
- Autoclavable.
- Colors in assorted packaging: Blue, green, red, silver, white, black, orange and pink.



REG. ANVISA: Dispensa

VARIOUS

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698761050	Hexagonal Instrument Marker - Packaging containing 120 units in Black - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698761043	Hexagonal Instrument Marker - Packaging containing 120 units in Orange color - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698761036	Hexagonal Instrument Marker - Packaging containing 120 units in Grey color - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698761029	Hexagonal Instrument Marker - Packaging containing 120 units in White color - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698761012	Hexagonal Instrument Marker - Packaging containing 120 units in Blue color - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698761005	Hexagonal Instrument Marker - Pack containing 120 units in Assorted color - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698760992	Hexagonal Instrument Marker - Pack containing 120 units in Pink - WITH ORGANIZING PACKAGING	10	-	3917.32.40
7898698760985	Hexagonal Instrument Marker - Pack containing 120 units in Red - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760978	Hexagonal Instrument Marker - Pack containing 60 units in Green - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760961	Hexagonal Instrument Marker - Pack containing 60 units in Black - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760954	Hexagonal Instrument Marker - Pack containing 60 units in Orange - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760947	Hexagonal Instrument Marker - Pack containing 60 units in Gray color - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760930	Hexagonal Instrument Marker - Pack containing 60 units in White - WITH ORGANIZING PACKAGING	10	-	340700.20

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760923	Hexagonal Instrument Marker - Packaging containing 60 units in Black - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760916	Hexagonal Instrument Marker - Packaging containing 60 units in Orange color - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760909	Hexagonal Instrument Marker - Packaging containing 60 units in Grey color - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760893	Hexagonal Instrument Marker - Packaging containing 60 units in White color - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760886	Hexagonal Instrument Marker - Packaging containing 60 units in Blue color - WITH ORGANIZING PACKAGING	10	-	340700.20
7898698760879	Hexagonal Instrument Marker - Pack containing 100 units in Assorted color	10	-	340700.20
7898698760862	Hexagonal Instrument Marker - Pack containing 100 units in Pink	10	-	340700.20
7898698760855	Hexagonal Instrument Marker - Pack containing 100 units in Red	10	-	340700.20
7898698760848	Hexagonal Instrument Marker - Pack containing 100 units in Green	10	-	340700.20
7898698760831	Hexagonal Instrument Marker - Pack containing 100 units in Black	10	-	340700.20
7898698760824	Hexagonal Instrument Marker - Pack containing 100 units in Orange	10	-	340700.20
7898698760817	Hexagonal Instrument Marker - Pack containing 100 units in Gray color	10	-	340700.20
7898698760800	Hexagonal Instrument Marker - Pack containing 100 units in White	10	-	340700.20

POTE DAPPEN

Glass Dappen Pot

RECOMMENDATION:

- Pot for handling materials.

FEATURES AND ADVANTAGES:

- It has two cavities for manipulation;
- Autoclavable 121 C;
- Colors: amber, blue, green and colorless.

REG. ANVISA: Dispensa



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468703	Amber Glass Dappen Pot	12	—	7017.20.00
7898936468727	Blue Color Glass Dappen Pot	12	—	7017.20.00
7898936468710	Green Color Glass Dappen Pot	12	—	7017.20.00
7898936468697	Colorless Dappen Glass Jar	12	—	7017.20.00

REG. ANVISA: Dispensa



POTE PALADON

Glass Paladon pot

RECOMMENDATION:

- Pot for handling materials.

FEATURES AND ADVANTAGES:

- With tampa;
- Autoclavable 121 C;
- Colorless.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468024	Paladon Pot	6	—	7017.90.00

REG. ANVISA: 80442020021

PULPOSUL

Dental cement

RECOMMENDATION:

• Sedative base, cavity liner, sealing dressing and temporary cementation of crowns, bridges, inlays and temporary fillings, when removals are planned, serve as an intermediate base under permanent restorations.

FEATURES AND ADVANTAGES:

- Resistant to friction;
- Good adhesion to the cavity walls;
- Poor solubility;
- Easy removal.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468031	Pulposul Liquid	12	—	3006.40.12
7898936468048	Pulposul Powder	12	—	3006.40.12
7898936469779	Pulposul kit	6	—	3006.40.12

REG. ANVISA: 80442020060

CIMENTO CIRÚRGICO

Dental cement

RECOMMENDATION:

Indicated as a periodontal coating to protect sutures from periodontal surgeries.

FEATURES AND ADVANTAGES:

- Excellent adhesion and plasticity properties;
- Safe local protection for surgical wounds, keeping them clean.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467812	Liquid Surgical Cement	12	—	3006.40.12
7898936467829	Surgical Cement Powder	12	—	3006.40.12
7898936469762	Surgical Cement kit	6	—	3006.40.12

VARIOUS

ALVEOLIPAST

REG. ANVISA: 80442029007

RECOMMENDATION:

- It is a paste recommended for treatments of alveolar infections (alveolitis) and helping with tissue regeneration.

FEATURES AND ADVANTAGES:

- Excellent healing agent;
- Free from Eugenol;
- Helps with post-surgical recovery,





BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936460727	Alveolipast 3g syringe	6	—	3006.40.12
7898936467737	Alveolipast 10g bottle	12	—	3006.40.12

WHITE DENT HOME

Tooth whitening

REG. ANVISA: 80442020079

RECOMMENDATION:

- Whitening agent based on 10%, 16% and 22% carbamide peroxide for home use, supervised by the dentist.

FEATURES AND ADVANTAGES:

- Neural action of potassium nitrate and occlusive action of sodium fluoride. Potassium nitrate is an effective agent in the treatment of hyperesthetic dentin. Its action is based on the fact that the increase in the concentration of extracellular potassium depolarizes the membranes of nerve fibers, blocking the passage of the stimulus, reducing hypersensitivity.
- Does not drip in the tray.






BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469465	1-WHITEDENT HOME 10% - Pack containing 1 syringe with 3g of product and 1 tip.	10	--	3306.90.00
7898936469472	2- WHITEDENT HOME 10% - Cartridge containing 3 syringes with 3g of product and 3 tips	10	--	3306.90.00
7898936469496	3 -WHITEDENT HOME 10% - Box containing 5 syringes with 3g of product, 5 tips, 2 trays and 1 tray storage case	10	--	3306.90.00
7898936469502	4 - WHITEDENT HOME 10% - Cartridge containing 10 syringes with 3g of product and 10 tips	10	--	3306.90.00
7898936469519	5 - WHITEDENT HOME 16% - Pack containing 1 syringe with 3g of product and 1 tip.	10	--	3306.90.00
7898936469526	6 - WHITEDENT HOME 16% - Cartridge containing 3 syringes with 3g of product and 3 tips	10	--	3306.90.00
7898936469540	7 - WHITEDENT HOME 16% - Box containing 5 syringes with 3g of product, 5 tips, 2 plates for trays and 1 case for storing trays	10	--	3306.90.00
7898936469557	8 - WHITEDENT HOME 16% - Cartridge containing 10 syringes with 3g of product and 10 tips.	10	--	3306.90.00
7898936469564	9 - WHITEDENT HOME 22% - Pack containing 1 syringe with 3g of product and 1 tip.	10	--	3306.90.00
7898936469571	10-WHITEDENT HOME 22% - Cartridge containing 3 syringes with 3g of product and 3 tips	10	--	3306.90.00

RELEASE

REG. ANVISA: 80442029014



DYCALCIUM

RECOMMENDATION:

- Recommended for pulp capping and protective lining under restorative materials, cements and other base. Biocompatible, does not contain eugenol does not inhibit setting of self-curing or light-curing resins. Has color similar to dentin. Easy to mix.

FEATURES AND ADVANTAGES:

- Safety for the procedure: important as pulp protection, especially in deep cavities.
- It can be used in direct or indirect pulp capping and in lining cavities with excellent levels of absence of postoperative pain.
- Does not contain eugenol.
- High compressive strength: greater tensile strength than temporary cements on the market.
- Fast setting: fast hardening.

DYCALCIUM

Calcium Hydroxide Lining Cement



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760411	Dycalcium	5	—	3006.40.11

REG. ANVISA: 80442020057

MASTER DAM

Light-cured Gingival Barrier



RECOMMENDATION:

- Used to protect the gingival tissue, in clinical procedures that may cause discomfort or irritation due to contact with the products used.

FEATURES AND ADVANTAGES:

- Excellent gum isolation in in-office teeth whitening procedures;
- Light-curing in 30 seconds;
- Easy removal after tooth whitening.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469359	Master Dam green	6	—	3006.40.12
7898936469342	Master Dam blue	6	—	3006.40.12

WHITENERS

SENSITIVITY

Desensitizing gel 0.2% and 2%

RECOMMENDATION:

- Low viscosity gel based on potassium nitrate and sodium fluoride, indicated to alleviate tooth sensitivity.

FEATURES AND ADVANTAGES:

- Transparent gel;
- Double desensitizing action, through potassium nitrate and sodium fluoride;
- Easy application.

REG. ANVISA: 80442020059



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467380	Sensitivity 0,2%	6	—	3006.40.12
7898936467397	Sensitivity 2%	6	—	3006.40.12

REG. ANVISA: 80442020036

MOLDEIRA TERMOMOLDÁVEL

Thermomoldable tray for whitening

RECOMMENDATION:

Thermo-moldable tray for whitening, making sports or orthodontic mouthguards.

FEATURES AND ADVANTAGES:

- Composed of ethylene vinyl acetate;
- No special equipment is required to make it;
- It is molded using only hot water.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469038	Tray (pair)	6	—	3926.90.40

REG. ANVISA: 80442020062



PLACA WHITEDENT

Whitening and bruxism plate

RECOMMENDATION:

Recommended for making wood for teeth whitening and bruxism.

FEATURES AND ADVANTAGES:

- Composed of ethylene vinyl acetate;
- Available in versions with square or round plate, compatible with any laminator;
- Biocompatible with the organism.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469984	Square 1mm - 2 pcs	12	—	3901.30.90
7898698760008	Square 1mm - 5 pcs	12	—	3901.30.90
7898698760046	Square 2mm - 2 pcs	12	—	3901.30.90
7898698760084	Square 3mm - 2 pcs	12	—	3901.30.90

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469991	Square 1mm - 2 pcs	12	—	3901.30.90
7898698760015	Square 1mm - 5 pcs	12	—	3901.30.90
7898698760053	Square 2mm - 2 pcs	12	—	3901.30.90
7898698760091	Square 3mm - 2 pcs	12	—	3901.30.90

Iodontosul always seeking the best for its customers and professionals who use its products sought to better understand their needs, evaluate suggestions and criticisms through "feedback" with renowned dental surgeons about the ICE TEST product.

Dra. Flávia Baldasso, CRO/RS 20369



"Like any good endodontist, I chose to work with just one brand, which was Iodontosul, I use their line of endodontic products because it is very complete, from sensitivity testing, to irrigating solutions and fillings. ICE TEST came fill a gap in the market, providing quality at an affordable price. I am a university professor in the city of Caxias do Sul and we know that the dentistry course is not cheap. My suggestion for students is to always purchase products that have quality and a good price. In this sense, Iodontosul is at the forefront. In addition to making it easier to purchase its products, ICE TEST has a quality that few have, when applying the product jet to the cotton, it does not spread into the environment, it is completely targeted. Furthermore, it preserves its temperature until we take the product to the patient's tooth. I use and recommend Iodontosul products.»



IODONTOSUL
PRODUCTS FOR HEALTH

ICE TEST

Vitality test spray

REG. ANVISA: 880442020032 - Ice Test

RECOMMENDATION:

Recommended for diagnosing pulp vitality and sensitivity.

FEATURES AND ADVANTAGES:

- Contains cooling agent (-50°)
- Minty aroma.



ENDODONTIC

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468253	Ice Test	12	-	2711.1990

ÓLEO DE LARANJA

Gutta percha solvent

REG. ANVISA: 80442029013

RECOMMENDATION:

Agent used to empty the root canal.

FEATURES AND ADVANTAGES:

- It is an essential oil,
- Extracted from sweet orange peel;
- With a pleasant odor;
- Volatile,



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760404	Orange oil	12		-7017.20.00

REG. ANVISA: 80442020023



EUCALIPTOL

Eucalyptus oil

RECOMMENDATION:

Solvent based on eucalyptus essential oil for removing gutta from the pulp chamber and root canal.

FEATURES AND ADVANTAGES:

- Antibacterial action;
- Insoluble in water;
- High degree of purity;
- Fast action.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467997	Eucaliptol	12		2932.99.11

REG. ANVISA:

Solução: 80442020013
Gel: 80442029008

E.D.T.A TRISSÓDICO
Chelating agent



RECOMMENDATION:

• And a chelating agent for cleaning (removing the smear layer, dental mud) and widening the root canal for the adhesion of filling materials.

FEATURES AND ADVANTAGES:

- Easy application;
- Ready to use;
- Gel with adequate viscosity.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467904	E.D.T.A Solution	12	—	3006.40.12
7898698760572	E.D.T.A gel 24%	6	—	3006.40.12
7898698760787	E.D.T.A 100mL	6	—	3006.40.12
7898698760770	E.D.T.A 500mL	6	—	3006.40.12

REG. ANVISA: 80442020024

TERGENSUL
Channel detergent

RECOMMENDATION:

Wetting solution indicated for root canal irrigation.

ADVANTAGE FEATURE

- Acts as a physical cleaning agent;
- Non-irritating.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468116	Tergensul	12	—	2828.90.11

LÍQUIDO DAKIN

Sodium hypochlorite 0.5%

REG. ANVISA: 80442020014

RECOMMENDATION:

Solution indicated for instrumentation and irrigation of root canals, acting as a disinfectant.

FEATURES AND ADVANTAGES:

- Disinfectant;
- Bactericidal;
- Fungicide;
- Low surface tension.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467478	Dakin 1000ml	18	—	2828.90.11

REG. ANVISA: 80442020015



SOLUÇÃO DE MILTON

Sodium hypochlorite 1%

RECOMMENDATION:

Solution indicated for instrumentation and irrigation of root canals, acting as a disinfectant.

FEATURES AND ADVANTAGES:

- Disinfectant;
- Bactericidal;
- Fungicide;
- Low surface tension.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468109	Milton 1000ml	18	—	2828.90.11

REG. ANVISA: 80442029012

SOLUÇÃO DE LABARRAQUE

Sodium hypochlorite 2.5%



RECOMMENDATION:

Solution indicated for instrument and root canal irrigation, acting as a disinfectant.

FEATURES AND ADVANTAGES:

- Disinfectant;
- Bactericidal;
- Fungicide;
- Low surface tension.

CÓD. DE BARRAS	DESCRIÇÃO	UN	I.P.I	NCM
7898936468987	Labarraque 1000ml	18	—	2828.90.11

ENDODONTIC

SODA CLORADA

Sodium hypochlorite 5%

REG. ANVISA: 80442020018



RECOMMENDATION:

Solution indicated for instrumentation and irrigation of root canals acting as a disinfectant.

FEATURES AND ADVANTAGES:

- Disinfectant;
- Bactericidal;
- Fungicide;
- Low surface tension.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468093	Chlorinated soda 1000ml	18	—	2828.90.11

ATTACK ÁCIDO

Hydrofluoric conditioning acid 5% and 10%

REG. ANVISA: 80442020061

5%

10%



RECOMMENDATION:

Attack acid 5% and 10% is a hydrofluoric acid-based gel for surface conditioning of dental porcelain, helping to promote retention in restorations.

FEATURES AND ADVANTAGES:

- Presence of dye, facilitating identification of the conditioned surface;
- Excellent viscosity, without risk of runoff;
- Excellent affinity with water, allowing easy removal.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469328	Attack Ácido 5%	6	—	3006.40.12
7898936469335	Attack Ácido 10%	6	—	3006.40.12

ATTACK ÁCIDO

Phosphoric conditioning acid 37%

RECOMMENDATION:

Phosphoric acid-based gel for surface conditioning of enamel and/or dentin, allowing greater and lesser retention of resinous materials.

FEATURES AND ADVANTAGES:

- Presence of blue dye, facilitating identification of the conditioned surface;
- Excellent viscosity, without risk of runoff;
- Excellent affinity with water, allowing easy removal.



37%

REG. ANVISA: 80442029003

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469311	Attack Ácido 37%	50	—	3006.40.12

REG. ANVISA: 80442020010

FLUORSUL PASTILHA

Fluoride tablets

RECOMMENDATION:

Fluoride tablets indicated for obtaining a mouthwash solution and preventing tooth decay and dental remineralization.

FEATURES AND ADVANTAGES:

- High degree of purity;
- Diversity to obtain formulations;
- Does not contain dye;
- Does not contain preservatives.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467294	Insert 60 units	12	—	2826.1990
7898936467287	Insert 120 units	12	—	2826.1990

REG. ANVISA: 80442020006



REPLASUL P

Bacterial plaque revealing tablet

RECOMMENDATION:

Basic fuchsin-based tablet, indicated for identifying bacterial plaque on teeth and adjacent structures.

FEATURES AND ADVANTAGES:

- 1 bottle with 60 units;
- 1 bottle with 120 units.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468789	Replasul 60 units	12	—	3204.100
7898936468789	Replasul 120 units	12	—	3204.100

REG. ANVISA: 80442020044

FIO DENTAL

Floss for dental cleaning

RECOMMENDATION:

Helps reduce gum problems, interproximal and subgingival cavities by removing plaque between teeth and below the gum line.

FEATURES AND ADVANTAGES:

Glides more easily between teeth without fraying.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468475	Mint 25 meters	12	S.T.	3306.20.00
7898936468529	Mint 50 meters	12	S.T.	3306.20.00
7898936468574	Mint 100 meters	12	S.T.	3306.20.00
7898936468512	Tangerine 25 meters	12	S.T.	3306.20.00
7898936468567	Tangerine 50 meters	12	S.T.	3306.20.00
7898936468611	Tangerine 100 meters	12	S.T.	3306.20.00
7898936468505	Strawberry 25 meters	12	S.T.	3306.20.00

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468550	Strawberry 50 meters	12	S.T.	3306.20.00
7898936468604	Strawberry 100 meters	12	S.T.	3306.20.00
7898936468482	Mint 25 meters	12	S.T.	3306.20.00
7898936468536	Mint 50 meters	12	S.T.	3306.20.00
7898936468581	Mint 100 meters	12	S.T.	3306.20.00
7898936468499	Grape 25 meters	12	S.T.	3306.20.00
7898936468543	Grape 50 meters	12	S.T.	3306.20.00
7898936468598	Grape 100 meters	12	S.T.	3306.20.00

PREVENTION

REPLASUL S

Bacterial plaque revealing solution

RECOMMENDATION:

Basic fuchsin-based solution, used to highlight the presence of bacterial plaque through mouthwash.

FEATURES AND ADVANTAGES:

Helps patients' oral hygiene;
Enables detection of bacterial plaques;
Interproximal.

REG. ANVISA: 80442020003



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468055	Replasil S	24	-	3204.13.00

CARIOSTASUL

Cariostatic

RECOMMENDATION:

Used to inhibit and prevent dental caries, especially in children, in the prevention of secondary cavities.

FEATURES AND ADVANTAGES:

- Silver nitrate at 10%, 12% and 30%;
- Strengthens the enamel structure while preventing the progression of cavities

REG. ANVISA:

80442020016 - 10%


80442020011 - 12%

80442020012 - 30%



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467782	Cariostasul 10%	12	-	3006.40.12
7898936467799	Cariostasul 12%	12	-	3006.40.12
7898936467805	Cariostasul 30%	12	-	3006.40.12

REG. ANVISA: 80442020082



MOLDEIRA DE CERA

Disposable Wax Tray

RECOMMENDATION:

Recommended for applying fluoride.

FEATURES AND ADVANTAGES:

Disposable mold made of malleable wax, covered with hydrophilic foam. Adaptable to the dental arch, facilitating the application of fluoride in the interproximal spaces.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698761128	Blue Wax Tray 24 units	10	-	340700.20
7898698761098	Blue Wax Tray 50 units	10	-	340700.20
7898698761098	Blue Wax Tray 100 units	10	-	340700.20
7898698761135	Yellow Wax Tray 24 units	10	-	340700.20
7898698761081	Yellow Wax Tray 50 units	10	-	340700.20

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760671	Yellow Wax Tray 100 units	10	-	340700.20
7898698761142	Assorted Wax Tray 24 units	10	-	340700.20
7898698761104	Assorted Wax Tray 50 units	10	-	340700.20
7898698761111	Assorted Wax Tray 100 units	10	-	340700.20

REG. ANVISA: 80442020005

REPLASUL T

Bacterial plaque revealing tablets

RECOMMENDATION:

Basic fuchsin-based solution, used to highlight the presence of bacterial plaque through topical application.

FEATURES AND ADVANTAGES:

- Helps patients' oral hygiene;
- Enables detection of bacterial plaque;
- Interproximal.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468062	Replasil T	12	-	3204.13.00

REG. ANVISA: 80442020062



REV CARIE

Cavity Viewer

RECOMMENDATION:

Solution used to highlight decayed tissue.

FEATURES AND ADVANTAGES:

Enables the selective removal of decayed tissue.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468161	Caries Viewer	12	-	3204.13.00

FLUORSUL GEL

Fluoride gel

RECOMMENDATION:

Gel indicated for preventing tooth decay and dental remineralization.
Acidified 1.23% and neutral 2%.

FEATURES AND ADVANTAGES:

- It has a great consistency for application;
- Decrease in the number and potential of oral microorganisms;
- Decreased tooth sensitivity;
- Promoting the remineralization of compromised tooth enamel;
- Aromas: mint, strawberry, tutti-frutti and unflavored.

REG. ANVISA: 80442020002 - Acidulado
80442020007 - Neutro

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467201	Acidulated Mint	12 / 50	—	2826.1990
7898936467218	Acidulated Strawberry	12 / 50	—	2826.1990
7898936467232	Acidulated Tutti Frutti	12 / 50	—	2826.1990
7898936467225	Acidulated without taste	12 / 50	—	2826.1990
7898936467249	Neutral Mint	12 / 50	—	2826.1990
7898936467256	Neutral Strawberry	12 / 50	—	2826.1990
7898936467270	Neutral Tutti Frutti	12 / 50	—	2826.1990
7898936467263	Neutral without taste	12 / 50	—	2826.1990



FLÚOR PLAX DAY

REG. ANVISA: 80442020069 Flúor plax Day

RECOMMENDATION:

Solution indicated for preventing tooth decay.

FEATURES AND ADVANTAGES:

- Decrease in the number and potential of oral microorganisms;
- Aromas: mint, strawberry, tutti-frutti and unflavored;
- The flavorless solution also has no dye, ideal for teeth whitening patients.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467126	Daily Flúorsul 1000mL - Mint	18	—	2826.1990
7898936467140	Daily Flúorsul 1000mL - Strawberry	18	—	2826.1990
7898936467164	Daily Flúorsul 1000mL - Unflavored	18	—	2826.1990
7898936467188	Daily Flúorsul 1000mL - Tutti Frutti	18	—	2826.1990

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467133	Plax day Fluoride 500mL - Mint	24	—	2826.1990
7898936467157	Fluoride Plax day 500mL - Strawberry	24	—	2826.1990
7898936467171	Plax day Fluoride 500mL - White	24	—	2826.1990
7898936467195	Fluoride Plax day 500mL - Tutti Frutti	24	—	2826.1990



FLÚOR PLAX WEEK

REG. ANVISA: 80442020067 Flúor plax Week

RECOMMENDATION:

0.2% Sodium Fluoride solution for weekly mouthwashes.

FEATURES AND ADVANTAGES:

Decrease in the number and potential of oral microorganisms.
Aromas: mint, strawberry, tutti-frutti and unflavored;
The flavorless solution also has no dye, ideal for teeth whitening patients.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467317	Weekly Flúorsul 1000mL - Mint	18	-	2826.1990
7898936467331	Weekly Flúorsul 1000mL - Strawberry	18	-	2826.1990
7898936467355	Flúorsul Weekly 1000mL - Unflavored	18	-	2826.1990
7898936467379	Weekly Flúorsul 1000mL - Tutti Frutti	18	-	2826.1990

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898698760312	Plax Week Fluoride 500mL - Mint	24	-	2826.1990
7898698760379	Plax Week Fluoride 500mL - Strawberry	24	-	2826.1990
7898698760336	Plax Week Fluoride 500mL - White	24	-	2826.1990
7898698760343	Plax Week Fluoride 500mL - Tutti Frutti	24	-	2826.1990



X RAY

REVEAL AND FIXER

REG. ANVISA: Dispensa

X Ray R Developer

RECOMMENDATION:

The product is indicated as an image developer in dental x-rays.

RADIOLOGY

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469427	fixative	12	-	370790.10
7898936469410	developer	12	-	370790.10



X Ray F

Fixative

RECOMMENDATION:

The best solution for fixing radiographic images on X-ray films.



ARMAÇÃO PLÁSTICA

Periapical x-ray storage

RECOMMENDATION:

Transparent plastic card for storing periapical radiographs.

FEATURES AND ADVANTAGES:

- Facilitates the viewing of x-rays, in addition to preserving them.
- Versions: 2 holes, 6 holes and 14 holes.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936467744	2-hole frame	10	—	9022.90.90
7898936467751	6-hole frame	10	—	9022.90.90
7898936467768	14-hole frame	10	—	9022.90.90

REG. ANVISA: 80442020014



ZNO

Zinc Oxide

REG. ANVISA: 80442020019

RECOMMENDATION:

Basic fuchsin-based solution, used to highlight the presence of bacterial plaque through topical application.

FEATURES AND ADVANTAGES:

- Helps patients' oral hygiene;
- Enables detection of bacterial plaque;
- Interproximal.



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468062	caries developer	12	—	3204.13.00

REG. ANVISA: 80442020017



EUGENOL

Temporary restorative material

RECOMMENDATION:

Solution used to highlight decayed tissue.

FEATURES AND ADVANTAGES:

Enables the selective removal of decayed tissue.

CÓD. DE BARRAS	DESCRIÇÃO	UN	I.P.I	NCM
7898936468161	Revelador de cárie	12	—	3204.13.00

MRI

Intermediate restorative material

REG. ANVISA: 80442020053

RECOMMENDATION:

Recommended for long-term temporary restorations.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468956	Liquid MRI 15ml	6	—	3006.40.11
7898936468949	MRI powder 38g	6	—	3006.40.11
7898936468932	MRI set	6	—	3006.40.11

HIGH PERFORMANCE:

- Temporary restorative cement;
- High dimensional stability;
- Low solubility;
- Durability.



CAVITEMP

Temporary shutter

REG. ANVISA: 80442029004 - Com Flúor
80442029006 - Livre de flúor

RECOMMENDATION:

Used in temporary cavity restoration procedures.

FEATURES AND ADVANTAGES:


- Provides adhesion to the cavity walls and a satisfactory seal;
- Prevents microleakage by organic fluids from the oral environment.



BAR CODE	DESCRIPTION	UN	I.P.I.	NCM
7898936469403	Cavitemp with fluoride	12	-	3006.40.11
7898936469366	Fluoride-free cavitemp	12	-	3006.40.11



FIBERGLASS PIN



Support Glass

Fiberglass pin

HIGH PERFORMANCE:

It has two formats: DOUBLE TAPER/RETAINER TAPER. Excellent adaptation to root canals, without the need for extensive wear of the tooth structure in the apical region.

REG. ANVISA: 8044209002




**DUPLA
CONICIDADE**



**SUPERFÍCIE
RETENTIVA**



BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468741	Support Glass nº 0,5 double conical	12	—	9021.2900
7898936467973	Support Glass nº 1 double conical	12	—	9021.2900
7898936468192	Support Glass nº 2 double conical	12	—	9021.2900
7898936468178	Support Glass nº 3 double conical	12	—	9021.2900
7898936468635	Support Glass nº 1 retainer conical	12	—	9021.2900
7898936468642	Support Glass nº 2 retainer conical	12	—	9021.2900
7898936468659	Support Glass nº 3 retainer conical	12	—	9021.2900



TARTARIT

Stain remover

RECOMMENDATION:

Recommended for removing stains caused by dental tartar.

FEATURES AND ADVANTAGES:

The product must be applied to the place where tartar and dental stains will be evident.

STAIN REMOVER



REG. ANVISA: 80442020055

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936468994	Tartarit	12	S.T	330690.00



HIGH PERFORMANCE

Lubricated parts can be sterilized without losing the lubricating and antioxidant properties of the oil.

BAR CODE	DESCRIPTION	UN	I.P.I	NCM
7898936469274	Dental Oil 100mL	12	ST	2710.1932
7898936467508	Dental Oil 200mL	12	ST	2710.1932
7898698760138	Triple applicator nozzle	5	-	8481.90.10

Dental Oil

Lubricating oil

REG. ANVISA: Dispensa

Triple applicator nozzle



IODONTOSUL
PRODUCTS FOR HEALTH

Us
follow
in the media
social

-  @iodontosul
-  @iodontosuloficial
-  www.iodontosul.com.br
-  +55 51 98594.2397

Rua Vitor Valpírio, 343 Bairro Anchieta
Porto Alegre/RS

The images in this catalog are for illustrative purposes only, and may change during the period of validity.